




# ECO 2



**Better Performance  
with Doppler Function**



# Why ECO?



# Key Benefits

- Affordable B&W system with PW function, extended diagnose way
- Advanced technologies allow **excellent image** quality;
- Ultra-compact & ultra-slim, easy to carry;
- Fashion design with **high resolution LED monitor**;
- Ultra-long battery life for long time use, chargeable from the car;
- **Easy to learn and use**: 6-one-key protocol , everything according to your intuition
- **Full Screen** function to enlarge the image;



- Hassle-free data management: all images and reports can be saved digital during the scan, and review on PC in a PC-readable format: jpg , PDF;
- Large storage:** 8G internal memory, data back up can be done in batch on external USB or PC effortless;
- Network:** press one key to send data directly to DICOM(as an option),DICOM worklist, storage ,print, structure report.
- Complete solutions for accessories



# Fashional Design



## ● 12 inch LED monitor

- Bright image even scan outdoor;
- Low energy consuming;
- Longer life than LCD;



## ● Elegant keyboard design

Ultra-compact  
Ultra-light  
Easy to carry







**Monitor is rotatable:**  
0 ~ 30°: tilt to different viewing angle

**Various output ports:**  
DICOM,VIDEO

**Dual probe holder**



# One solution with cart (TR9000)

## All-in-One

- Fashional design
- Powerful Function
- Unit can be fixed on the table through the screw
- (fasten the unit & theft-proof)



# Excellent Performance

- Multi-beam Former
- Space Compound Imaging
- SRA(speckle reduction technology )
- i-Image
- THI
- Chroma





# Intelligent Workflow

- Intelligent control panel design
- **Full screen Display**
- EasyView™
- Several one key function
- Full range calculation package
- User defined preset



# Intelligent Control Panel Design



- Independent areas for specific functions.
- Backlit keyboard can indicate active keys clearly.

Function key area




Printing area



Measurement menu



One key optimization



# Full Screen Display




**One key to enlarge the image to full screen , easy to view from distance.**



# Excellent PW

Excellent PW resolution and automatic PW trace measurement will greatly improve the diagnosis confidence by delivering both qualitative & quantitative results.



# EasyView™







# Language Support

- Chinese
- English
- German
- Spanish
- Russian
- Turkish
- French
- Czech
- Italian
- Portuguese




Change the language without restarting



# One key upgrade—easy to maintain



# Complete solution for accessories





Item	Part number	Picture
<b>Printer</b>	Sony 895MD Sony 897MD Sony 898MD Sony UPD711MD MISUBISHI P-93W	
	HP 1102 HP M251 HP 2055D HP 252dw HP 6230 EPSON L130	
<b>Dedicated Trolley</b>	TR9000	



Item	Part number	Picture
Carrying case	BG-100	
Battery	BT-2500	
Anti-water keyboard cover		
Vehicle Power Charger		



## Biopsy kits Configure

Part Mode	Probe type	Picture
BC04-A	C3-A	
BL03-A	L7M-A	
BT03-A	V6-A	
BC06-A	MC3-A	







**Thank you!**  
CHISON

We reserve the rights to make changes  
to this presentation without prior notice.

